

## Exercise no. 7

### Determination of water hardness

#### Determination of $\text{Ca}^{2+}$ and indirect determination of $\text{Mg}^{2+}$

##### *Introduction*

The ions involved in **water hardness**, i.e.  $\text{Ca}^{2+}$  and  $\text{Mg}^{2+}$ , can be determined by titration with a chelating agent, ethylenediaminetetraacetic acid (EDTA), usually in the form of disodium salt ( $\text{H}_2\text{Y}^{2-}$ ). Eriochrome Black T (Ind) is commonly used as indicator for the above titration. At pH 10,  $\text{Ca}^{2+}$  and  $\text{Mg}^{2+}$  ions first complex with the indicator as Ca-Ind and Mg-Ind which is **wine red**. As the stronger ligand EDTA is added, Ind from complexes is replaced by the EDTA, the colour is **blue**. The endpoint of titration is indicated by a sharp colour change from wine red to blue.

**PART A:** Titration using **Eriochrome Black T** (buffer solution) as indicator determines **total hardness due to  $\text{Ca}^{2+}$  and  $\text{Mg}^{2+}$  ions**.

**PART B:** **Hardness due to  $\text{Ca}^{2+}$  ion** is determined by a separate titration at a higher pH, by adding NaOH solution to precipitate  $\text{Mg}(\text{OH})_2\downarrow$ , using **murexide** as indicator.

##### *Experimental Procedures*

Dilute the analytical sample with distilled water in a measuring flask to the mark (100 mL) and mix completely. The solution in the measuring flask should be treated as 100 mL of the water sample.

### **Part A: Determination of total hardness**

**Pipette** 10 mL of the solution into the conical flask and add 5 mL of the ammonia buffer solution (pH=10). Swirl to mix. Do not add deionized water to the solution! Add one small pinch of Eriochrome Black T indicator, immediately prior to titrating a sample. Swirl to mix. The solution should be a pale pink. **Do not add more indicator to make the solution darker as this can cause problems with the endpoint.** Titrate the solution immediately with EDTA (0.01M) against a white background until the **light wine red** solution turns a **light sky blue**. Read the final volume (save the solution for colour comparison). Repeat the titration to obtain three concordant results (not differing more than 0.2 mL).

### **Part B: Determination of concentration of Ca<sup>2+</sup> ions**

**Pipette** 10 mL of the water sample (from 100 mL) to the Erlenmeyer flask and add 10 mL of 1 M NaOH. Swirl to mix. Do not add deionized water to the solution! Add one small pinch of **murexide** indicator, immediately prior to titrating a sample. Swirl the solution and wait for a one minute to completely precipitate the magnesium ions as Mg(OH)<sub>2</sub>↓. The solution should be a pale pink-red. **Do not add more indicator to make the solution darker as this can cause problems with the endpoint.** Titrate the solution immediately with EDTA (0.01M) against a white background until the **light pink-red** solution turns a **light blue-violet**. Read the final volume (save the solution for colour comparison). Repeat the titration to obtain three concordant results (not differing more than 0.2 mL).

### **Calculation**

1. From the results in **Part A**, calculate the hardness in °dH (German hardness unit, equivalent to CaO):

$$\text{°dH} = v_1 \cdot c_{\text{EDTA}} \cdot 56.08/10 \cdot 10$$

$v_1$  – volume of EDTA used in the presence of Eriochrome Black T indicator [mL]

2. From the results in **Part B**, calculate the concentration of  $\text{Ca}^{2+}$  ions in the water sample in mg/100mL:

$$m_{\text{Ca}} = v_2 \cdot c_{\text{EDTA}} \cdot 0.04008 \cdot 10$$

$v_2$  – volume of EDTA used in the presence of murexide indicator [mL]

$c_{\text{EDTA}}$  – concentration of EDTA [mmol/mL]

0.040008 – mass of a millimole of  $\text{Ca}^{2+}$  [g/mmol]

3. Hence, calculate the concentration of  $\text{Mg}^{2+}$  ions in the water sample in mg/100mL:

$$m_{\text{Mg}} = (v_1 - v_2) \cdot c_{\text{EDTA}} \cdot 0.02432 \cdot 10$$

$v_1$  – volume of EDTA used in the presence of Eriochrome Black T indicator [mL]

$v_2$  – volume of EDTA used in the presence of murexide indicator [mL]

$c_{\text{EDTA}}$  – concentration of EDTA [mmol/mL]

0.02432 – mass of a millimole of  $\text{Mg}^{2+}$  [g/mmol]